



Thermo Fisher Scientific

Isodat 3.0

Simultaneous ASCII Data Export

Revision A - 1250500



© 2008 Thermo Fisher Scientific Inc. All rights reserved.

All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Thermo Fisher Scientific Inc. provides this document to its customers with a product purchase to use in the product operation. This document is copyright protected and any reproduction of the whole or any part of this document is strictly prohibited, except with the written authorization of Thermo Fisher Scientific Inc.

The content of this document are subject to change without notice. All technical information in this document is for reference purposes only. System configurations and specifications in this document supersede all previous information received by the purchaser.

Thermo Fisher Scientific Inc. makes no representations that this document is complete, accurate or error-free and assumes no responsibility and will not be liable for any errors, omissions, damage or loss that might result from any use of this document, even if the information in the document is followed properly.

This document is not part of any sales contract between Thermo Fisher Scientific Inc. and a purchaser. This document shall in no way govern or modify any Terms and Conditions of Sale, which Terms and Conditions of Sale shall govern all conflicting information between the two documents.

Printing History: Revision A printed in November 2008.

Contents

Isodat 3.0 - Simultaneous ASCII Data Export.....	1
Activating Script Bar	2
Enabling Snapshot Script	3
Accessing Result Files	4
Structure of a Single Result File	6

Figures

Marking Script Bar checkbox	2
Selecting snapshot script	3
Close all open chromatograms	3
Selecting Results tab of File Browser	4
Contents of folder ACQ-Results	4
Contents of subfolder SnapShot	5
Contents of one subfolder Acquisition-000k	5
Structure of a result file	6
Close all open chromatograms	7

Isodat 3.0 - Simultaneous ASCII Data Export

This manual describes how raw data are exported by Isodat 3.0 as ASCII data (in csv format) during any running Continuous Flow acquisition. Depending on the acquisition time of the individual acquisition, directories will be created for saving fractions of the exported data.

In addition to the common separate data export after a measurement in Isodat 3.0, this gives the user online access to results already before an acquisition is finished. Thereby, data interpretation and measures for process control depending on signal intensity are possible promptly.


This manual deals with the following topics:

- [“Activating Script Bar” on page 2-2](#)
- [“Enabling Snapshot Script” on page 2-3](#)
- [“Accessing Result Files” on page 2-4](#)
- [“Structure of a Single Result File” on page 2-6](#)


Activating Script Bar



Start Isodat Acquisition.



Right-click on the Isodat Acquisition title bar.



Choose **Properties**.




Figure 1. Marking Script Bar checkbox

Mark the **Script Bar** checkbox. See [Figure 1](#).



The script bar will appear. Move the script bar until it is arranged in a row with the other bars.

Enabling Snapshot Script



Click on the **New** button.





Figure 2. Selecting snapshot script

Browse to the directory \Thermo\Isodat NT\Global\ISL\Snapshot. At Files of Type select **ISL Script**. See [Figure 2](#).

Double-click on the snapshot script **SwitchSnapShot.isl**. A new button, SwitchSnapShot, appears at the script bar.



To enable the snapshot script, click on the **SwitchSnapShot** button.



Confirm with **OK**.




Figure 3. Close all open chromatograms

Close all open chromatograms, especially the measurement file acquired at least. See [Figure 3](#).

Accessing Result Files

1. Start a sequence.
2. In the File Browser, select the **Results** tab.

The results folder ACQ-Results will be shown. See [Figure 4](#).





Figure 4. Selecting Results tab of File Browser

3. Right-click into the grid of the File Browser and select **Launch Explorer**.



As usual, the folder ACQ-Results contains chromatograms as dxf files. Their number equals the number of sequence lines you previously specified during sequence creation. Additionally however, a subfolder SnapShot has just been created. See [Figure 5](#).




Figure 5. Contents of folder ACQ-Results

For any of the acquisitions, a corresponding subfolder with the same name will be created in SnapShot.




Figure 6. Contents of subfolder SnapShot

See [Figure 6](#) as an example with three folders (Acquisition-0001, Acquisition-0002 and Acquisition-0003).




Figure 7. Contents of one subfolder Acquisition-000k

Each of these folders contains measured result files in csv format, that is as comma-separated values. The csv files are named 1.csv, 2.csv, etc. Long measurements will lead to more csv files than short ones. See [Figure 7](#).

Within each acquisition folder numbering starts anew with 1.csv, 2.csv, etc.. The individual export files (*.csv) can now be examined one after another in another program, for example a spreadsheet, while the measurement is still running.

Structure of a Single Result File

All csv-files have the same structure. Figure 8 shows it for an arbitrary csv file.

	A	B	C	D	E
1	200	115.541	134.399	148.128	
2	400	115.547	134.374	148.154	
3	600	115.541	134.373	148.154	
4	800	115.541	134.352	148.08	
5	1000	115.54	134.354	148.044	
6	1200	115.544	134.326	148.002	
7	1400	115.535	134.307	148.087	
8	1600	115.539	134.387	148.081	
9					
10	acquisition	intensity of each			
11	time [ms]	individual channel [mV]			
12					


Figure 8. Structure of a result file

The leftmost column A lists the acquisition times in ms. They are given in steps of ascending order which reflect the selected integration time.

Each of the following n columns (that is B, C, D,...) represents the intensities one of the n measured channels. As an example, three measured channels (m/z 44, m/z 45 and m/z 46) will lead to three columns.

Usually, a csv file comprises eight lines. In seldom cases however, it may contain more or less than eight lines. According to the default settings previously adjusted in Isodat Configurator, each intensity value is given with three decimal places.

Note Isodat will not delete any of the generated csv files. Delete the csv files you do not need anymore on your own. ▲



To disable the snapshot script, click on the **SwitchSnapShot** button.



Confirm with **OK**.



Figure 9. Close all open chromatograms

Close all open chromatograms. See [Figure 9](#). Creation of further csv files will be stopped.

Thermo Fisher Scientific Inc.

81 Wyman Street

P.O. Box 9046

Waltham, Massachusetts 02454-9046

United States

www.thermo.com